

# DRE Citadel EZ

## Electrosurgical Unit

### Ideal for light to moderate surgical procedures

The DRE Citadel EZ is a high frequency electrosurgical unit designed for light to moderate surgical procedures. Available with a wide range of electrosurgical accessories.

#### Features

- ✓ High frequency, 120-watt electrosurgical unit ideal for light to moderate surgical procedures.
- ✓ Multi-functional ESU operates in the following modes: Pure Cut, Blend, Forced Coag, Soft Coag, and Bipolar.
- ✓ Neutral electrode connection is constantly monitored for safe operation.
- ✓ Control the monopolar handpiece output for complete attention to the surgical field.
- ✓ Self check controls help to ensure all parts and accessories are functioning properly.
- ✓ Self patient plate circuit monitoring for patient safety.
- ✓ Self impedance monitoring for continuous output during surgical procedures.
- ✓ Also available in 160-watt configuration. Ask your DRE representative for details.

#### Standard Accessories

- ✓ Power cord
- ✓ Reusable electrode cable
- ✓ Reusable electrode plate
- ✓ Electrosurgical tips (pack of 10)
- ✓ Reusable hand switching pen
- ✓ Reusable foot switch
- ✓ Bipolar adapter



#### SPECIFICATIONS



**Weight:**  
11 lbs



**Dimensions:**  
10" x 4.1" x 11.3"



**Working Frequency:**  
600 KHz

**Selectable Input Voltage:**  
115 - 230 Vac

**Mains Frequency:**  
50 - 60 Hz

**Electrical Input Power:**  
300 VA



**CUT - max output:**  
120 W-250  $\Omega$

**BLEND - max output:**  
90 W-200  $\Omega$

**FORCED COAG - max output:**  
80W-150 $\Omega$

**SOFT COAG - max output:**  
60 W-100  $\Omega$

**BIPOLAR - max output:**  
40 W-100  $\Omega$

1547-01-2017-06