

SERVING CUSTOMERS IN  
**OVER 150 COUNTRIES**

Supplying  
new and  
refurbished  
medical  
equipment  
since 1984



- » **Anesthesia**
- » **Surgical Tables**
- » **Lighting**
- » **Stretchers**
- » **Monitoring**
- » **Ventilators**
- » **Infant Care**
- » **Defibrillators**
- » **Electrosurgical**

*and much more  
for the entire O.R.*

Introducing the  
**DRE Ventura CS2**  
Anesthesia System  
Page 4



# Operating within your budget for more than 25 years



**DRE has been a trusted leader in the international medical equipment industry since 1984.**

We provide new and professionally refurbished medical equipment to clients in more than 150 countries. Our extensive customer base includes hospitals, surgery centers, mission organizations, and distributors.

**One of many ways that DRE can help you is by serving as the single source for virtually all of your medical equipment needs.**

Our valued clients save money and time by purchasing all of their equipment from DRE instead of paying multiple independent vendors for both products and shipping expenses. Our 18,000 square foot facility is located in Louisville, Kentucky (USA) and is stocked with new and refurbished products that can be shipped to nearly every country. We can provide you with a wide range of products to furnish your operating theatre, intensive care unit, and patient recovery areas.

**Please call us at 502-244-6333 to learn more about how we can help you save time and money on your medical equipment purchases.** Our international

sales division is comprised of eleven experienced, multilingual professionals who have the knowledge to help you make the best purchasing decisions. We look forward to working with you to provide quality equipment solutions at the best possible price.




### Meet the DRE International team

Front row, left to right: Renata Wilkins, Amanda Cannady, Jasen Hargrove, Kary Van Arsdale, Sherri Evanoff, Kristin Easton. Back row, left to right: Barbara Young, Jay Thurman, Al Marquez, David Hardison, Francisco Cano

## See inside for new and refurbished medical equipment

Anesthesia .....	4
Stretchers .....	7
Surgery Tables .....	8
Lights .....	10
Monitors .....	12
Telemetry .....	13
Laparoscopy .....	13
Electrosurgical .....	14
Defibrillators .....	16
ECG/EKG Monitors .....	17
Crash Carts .....	17
Infant Care .....	18
Microscopes .....	20
Infusion Pumps .....	21
Suction Pumps .....	21
Ventilators .....	22




See our complete selection at  
**www.dremed.com**

- Browse thousands of new and professionally refurbished medical and surgical products
- Quickly search by equipment type, keyword, manufacturer and Web ID number
- Place a quote request for the equipment you need. We can ship to nearly any destination in the world!



**DRE Ventura CS2**

Modern anesthesia delivery at an affordable price



- High quality U.S. and European technology and components
- Advanced ventilation: PCV, PSV, SIMV, SMMV, volume control, and auto-peep
- Adult, pediatric, and neonatal modes
- Integrated digital oxygen monitor
- Three gas machine
- Autoclavable absorber (excluding manometer)

« **DRE Integra AV-S**

Advanced, user-friendly anesthesia system



- Standard with three gases (O<sub>2</sub>, N<sub>2</sub>O, air)
- Air/N<sub>2</sub>O interlock
- Open architecture gives you ample space for monitors and accessories
- Holds two vaporizers

**DRE Integra SP Series** »

Affordable and low maintenance

- Choose a system that holds one or two vaporizers
- Standard with three gases (O<sub>2</sub>, N<sub>2</sub>O, and Air)
- Variety of ventilator options available
- Air/N<sub>2</sub>O Interlock




**MRI-COMPATIBLE** configuration available

Compact, portable system that is great for mobile anesthesiologists

« **DRE Integra SP VS02**

- Three gas system (O<sub>2</sub>, N<sub>2</sub>O and air)
- Holds one vaporizer (your choice)
- Weighs only 29 pounds without vaporizer
- Designed for low flow anesthesia
- Anti-hypoxic guard
- Air/N<sub>2</sub>O interlock



Space saving design

Available with optional carrying case for safe transport

Also available with optional mobile cart

**Datex-Ohmeda Aestiva/5** »

Provides low and minimal flow anesthesia delivery

- Includes 7100 or 7900 series ventilator
- Innovative patient breathing system
- Automatic fresh gas flow compensation

Also available: **Datex-Ohmeda S/5 Aespire**



Web ID: 1469



Web ID: 59

Certified Refurbished

« **Datex-Ohmeda Excel 210 SE**

Ventilation and monitoring at an affordable price

- Features three gases and two vaporizers
- Integrated oxygen, volume and airway pressure monitoring

Also available: **Excel 110 and 210**



Web ID: 1670

« **Dräger Fabius GS**

Reliable, advanced ventilation components

- Electronically controlled E-vent
- Promotes easy monitoring of gas flow
- Intelligent care management

Also available: **Dräger Narkomed GS**

Certified Refurbished

**Dräger Narkomed GS** »

Excellent quality for facilities with limited budgets

- Three gas machine with integrated monitoring system
- Holds up to three vaporizers
- Includes the dependable AV2+ ventilator

Also available: **Dräger Narkomed 2B, 2C, 3 and 4**



Web ID: 1936

**DRE Millennium 5** »

Emergency stretcher with a fifth wheel for improved steering

- Back rest and knee rest lifting
- Central braking system
- 500-pound weight capacity
- High density mattress
- Two IV poles



Web ID: 2089

Easy to position stretcher with back and knee-rest lifting, and Trendelenburg



**Positioning**

- Head elevation: 0 - 75°
- Knee elevation: 0 - 30°
- Height adjustment: 573 mm - 925 mm

Well-placed handles and a fifth wheel for improved steering through your facility

Visit [www.dremed.com](http://www.dremed.com) for our complete selection of remanufactured stretchers and hospital beds!



**Hill-Rom 894 Century CC**

Web ID: 1857



**Hill-Rom Centra Series 850/852**

Web ID: 1860



**Hill-Rom Total Care**

Web ID: 1862

Serves as a  
Surgery Table,  
Procedure Chair,  
Stretcher and  
Recovery bed



Web ID:  
2183

**DRE Versailles P100**

- Power operated with hand remote control
- Full C-arm compatibility
- Traversing patient platform
- Lateral tilt with excellent stability
- Long life battery backup
- 550-pound weight capacity
- Includes deluxe mattress and IV pole



Variety of optional  
headrests available

Transport

Procedures

Surgery

Recovery



Ophthalmology



Gynecology



C-arm Coverage

Also available in a cost-effective manual configuration

Web ID:  
2199



**DRE Milano E20**

Power Procedure Chair

Web ID:  
2200



**DRE Milano T50**

Power Procedure Table

Web ID:  
2201



**DRE Verona**

Power Procedure Chair

**DRE Lucerne ES**

Versatile  
table provides  
advanced  
positioning

- Electro-hydraulic operation with remote control
- One button repositioning
- Horizontal slide feature
- Great for X-ray and C-arm use
- Optional orthopedic package available



Web ID:  
2194

**DRE Marseille**

Dependable, cost-effective  
manual table

- Foot pedal operation
- Trendelenburg and Reverse Trendelenburg positioning
- Easy to operate and maintain



Web ID:  
1669

**Amsco 3085**

Durable surgery table supports  
up to 1,000 pounds

- Electric operation with remote control
- Exceptional flexibility and weight capacity
- Built in safety and reliability

Also available  
Amsco 3080  
Supports up  
to 500 lbs.

Web ID:  
1499



Certified Refurbished

**Amsco 2080**

Affordable, dependable surgery table

- Available in manual, electric, and remote control versions
- Positions include: Raise/lower, lateral tilt, Trendelenburg, flex
- 300-pound weight capacity



Web ID:  
164

**DRE Vision LED** »


Efficient surgery light built to last for years



**Bright Light,  
Long Life,  
Low Heat**

Web ID: 2119

- Available in a variety of ceiling mount and portable styles
- 65,000 lux per light head
- LED technology lasts up to 50,000 hours
- 5-stage dimming
- Half the electric consumption of halogen lights



**DRE Maxx Luxx II** »

- Available in single, dual, and triple head styles
- Choose from ceiling mount, wall mount, and mobile styles
- 102,000 lux per light head
- Automatic bulb change mechanism



Web ID: 128


**Cool lighting system available with a video camera**



**Also available**



**DRE Maxx Luxx 3**  
Efficient Surgery Lights




**DRE Maxx Luxx Select**  
137,000 lux per light head


« **DRE Vision Excel**

Surgery light provides improved shadow control

- Single, dual, and mobile versions
- 64,500 lux per light head
- Large, pre-focused pattern
- 16" diameter elliptical reflector



Web ID: 1252



Vision Excel (Dual)



Vortex/Vision Excel



Vista/Vision Excel

« **DRE Vista**

Exam Light

- 43,000 lux
- Available in single, dual, and mobile versions
- Perfect for office examinations



Web ID: 1254

**DRE Vista LED** »

LED Exam Light

- Reduced energy consumption
- Long life expectancy



Web ID: 2196

Visit [www.dremed.com](http://www.dremed.com) for our complete selection of new and refurbished surgery lights!



**Skytron Infinity**

Web ID: 1804



**Amsco Quantum SQ-240**

Web ID: 125



**Amsco Polaris**

Web ID: 124



« **DRE Waveline Pro** Web ID: 1821

**Anesthesia Monitor with up to 8 waveforms**

- 15" high resolution, flat panel display
- Variety of configurations available
- Basic models include NIBP, SpO<sub>2</sub> and EKG
- Optional: CO<sub>2</sub>, IBP, Temp. and Agent Analysis



» **DRE Waveline Touch**

Web ID: 1979

**User-friendly touch-screen patient monitor**

- NIBP, SpO<sub>2</sub>, EKG and Temperature
- Also available with IBP and/or CO<sub>2</sub>
- Optional gas analysis module also available



Web ID: 1980

« **DRE Waveline EZ**

**Durable transport monitor with touch-screen**

- NIBP, SpO<sub>2</sub>, EGK and Temperature
- Also available with CO<sub>2</sub>
- Long-life battery (up to 8 hours!)

Visit [www.dremed.com](http://www.dremed.com) for our complete selection of new and refurbished patient monitors!



**Datex-Ohmeda S/5**

Web ID: 1706



**Datex-Ohmeda CardiCap/5**

Web ID: 1691



**Datascop Passport 2**

Web ID: 1690


» **DRE Envoy**

**Comprehensive central monitoring system**

- Monitors vital signs of multiple patients simultaneously
- Large, color touch-screen facilitates use of system
- Monitors from a single location with wireless technology
- Automatically compiles and archives discharge reports
- Transmits reports to hospital management systems



**Made in the USA**



**Easy-to-use interface**  
Review, edit and enter patient information while you observe all patients in the system.

**Laparoscopy & Endoscopy**

» **Stryker Laparoscopy Towers**

Certified Refurbished

**Customizable, all-in-one units**

- Choose from standard or HD models
- Camera options: Stryker 988, 1088 and 1188
- Light source options: Stryker X-6000, X-7000 and X-8000
- Insufflator options: Stryker 16L, 30L and 40L
- Stryker SCD Pro digital image capturer
- Choose from a variety of scopes and monitors



Web ID: 1866

See the difference **DIGITAL TECHNOLOGY** can make with **DRE ASG Series** electro-surgical units



**Features of the ASG Series include the following:**

- Cut, bipolar, coag, pinpoint, spray, and fulguration
- Ability to operate in either 110V or 220V areas
- Defibrillator-proof type CF
- Bright color-coded controls
- No calibration required
- Neutral electrode monitoring
- Optional carts available


**A VARIETY OF ACCESSORIES AVAILABLE**



Footswitch



Pen



Forceps



Grounding Cord/Plate



Web ID: 260

**DRE ASG 120**

**Ideal endoscopy unit**  
120 watts of cutting power.  
Two modes of coagulation.



Web ID: 258

**DRE ASG 200**

**For multi-purpose use**  
200 watts of cutting power.  
10 different blend settings.



Web ID: 244

**DRE ASG 300**

**For demanding procedures**  
300 watts of cutting power.  
Isolated RF output.

**Bovie IDS-400** >>

**400 watts of cutting power**

- Digital error detection
- Variety of advanced settings



Valleylab electro-surgical units for bipolar and monopolar procedures

**Valleylab Force FX**



Web ID: 421

- 300 watt output for cutting
- Three bipolar modes
- Instant response technology
- Automatically adjusts to skin tissue changes
- Improved safety and dependability



**OPTIONS**

- Cart (shown above)
- Foot pedal (Monopolar & Bipolar)
- Hand-switching pen
- Grounding pads (Disposable & Reusable)

Certified Refurbished

<< **Valleylab Force 2**



Web ID: 107

- Excellent combination of value and quality
- 300 watt output for cutting
- Standard adaptive REM system
- Permits 2 surgeons to operate simultaneously

Certified Refurbished

<< **Valleylab Force EZ**



Web ID: 1813

- Excellent performance at lower settings
- 300 watt output for cutting
- Instant response technology
- Two primary coagulation modes

Certified Refurbished

**Valleylab Ligasure** >>

- Effectively promotes vessel occlusion
- Fuses vessel walls to create permanent seal
- Available with a variety of pedals and instruments

Web ID: 1683



**Zoll M Series & R Series** »

Compact, portable, and bright

- Available new and refurbished
- Lightweight and portable
- Available with biphasic technology
- Optional parameters: NIBP, SpO<sub>2</sub>, AED and CO<sub>2</sub>



Web ID: 1293

« **Philips HeartStart XL**

Dependable, biphasic defibrillator

- Manual and Automatic (AED) operation
- Charges quickly
- Synchronized cardioversion

Web ID: 1487



« **HP Codemaster XL+**

Reliable general defibrillator

- Integrated EKG monitor and printer
- Built-in pediatric paddles
- Optional pacing

Certified Refurbished



Web ID: 90

« **Physio-Control Lifepak 12**

Portable defibrillator with options

- Reliable choice for ambulance or field use
- Biphasic or monophasic technology
- Optional parameters: NIBP, SpO<sub>2</sub>, AED and CO<sub>2</sub>

Web ID: 1435

Certified Refurbished

» **Physio-Control Lifepak 9**

Affordable basic defibrillator

- Integrated EKG monitor and printer
- Optional pacing
- 4" monitor display

Web ID: 87



« **DRE Velocity EKG**

- 12-channel EKG system
- Displays real-time EKG waveforms
- Large data storage
- Auto/manual measurement
- Includes analysis software
- Lightweight with carrying case



Web ID: 1947

Portable, PC-Based EKG System

**INSTALL**  
the software  
**CONNECT**  
the Velocity EKG  
**MEASURE**  
EKGs on your PC

Visit [dremed.com](http://dremed.com) for complete hardware requirements.




**DRE True ECG Advance**

12-channel EKG and stress test system

- Alpha-numeric keyboard
- One touch operation
- Options: Cart and treadmill

Web ID: 1977



**GE Mac 5500**

12-lead and 15-lead ECG systems

- For stress testing and resting EKG measurement
- Complete analytic algorithms
- Available with optional cart/stand

Web ID: 2103



Web ID: 1975

**DRE True ECG-3**

Three channels



Web ID: 1948

**DRE True ECG Plus**

Six channels



Web ID: 1988

**DRE True ECG 12**

Twelve channels

**DRE Crash Cart Packages**

Prepared by DRE to meet your emergency needs

- Basic units include IV pole and tank holder
- Ambu bags and suction pumps available
- 5 drawers with locking mechanism





« **DRE Infantia NB1**

**Infant incubator with servo control**

- Includes mattress, sensor, and IV pole
- LED digital display
- Touch switch panel
- Adjustable inclination of the bassinet
- Audible and visual alarms
- Set temperature, real temperature, and auxiliary temperature



« **Ohmeda Care Plus**

**Safe, secure infant incubator**

- Controls environment in patient and air modes
- Includes new patient mattress and sensor
- Multiple audio tones
- Instant response prevents unsteady temperature
- Multi-faceted alarm system (patient temperature, control temperature, high temperature, probe failure, air circulation, system failure and power failure)

Certified Refurbished

**DRE NeoGlider** »

**User-friendly transport incubator**

- Includes durable, collapsible mobile base
- Mattress and IV pole also included
- Double wall hood with side door
- Uses both AC power and DC power
- Displays set temperature, air temperature and baby temperature



**Air-Shields C-100** »

**Affordable incubator with digital display**

- Forced air circulation system
- Uniform heat distribution
- Digital temperature display
- Humidification
- Includes new patient mattress and sensor



« **DRE Neolite**

**Infant warming system available with optional phototherapy**

- Pre-warm, manual, and baby mode
- Adjustable infant bed
- Includes mattress, IV pole, and sensors
- Audible and visual alarms
- Displays set temperature and skin temperature



Optional phototherapy via LED lights

« **Air-Shields Isolette C550 QT-XL**

**Easy-to-use incubator with clear display**

- Controls for baby temperature and air temperature
- Includes new patient mattress and sensor
- Variable height adjustment
- Mattress tray tilt mechanism
- Provides stable temperature





Web ID: 2144

« **DRE Em5000**

**Fiber optic ENT microscope with premium European optics**

- Five step manual magnification changer
- Fiber optic illumination with cable
- 150 watt halogen light source
- Fine focus objective lens
- Wide-field binoculars
- Adjustable maneuvering handles
- Compact base with large casters



Web ID: 1464

**Includes a foot pedal for hands-free focus and zoom**

**DRE Om2100** »

**Automatic microscope perfect for ophthalmic and neurological procedures**

- Motorized focus and zoom
- 65 mm interchangeable objective lens
- Includes light source, fiber optic cable, and power supply
- Observation kit and beam splitter also included
- Sealed foot pedal
- Inclined binoculars
- Sturdy mobile stand with four-wheel base



Web ID: 184

« **Zeiss Opmi 6**

- Ideal for ophthalmic surgery
- Motorized focus and zoom
- Inclined binoculars with eyepieces
- High resolution
- Available with foot control

**Certified Refurbished**

**Zeiss Opmi 1 ENT** »

- Designed for ENT surgery
- Fine focus adjustment
- Straight binoculars
- Inclined coupling
- Light source with cable
- Heavy duty mobile stand

Web ID: 183



« **DRE Avanti M3**

Web ID: 1466

**Versatile pump accepts most types of straight line tubing**

- Lightweight and portable
- Perfect for use in hospitals or ambulances
- User-selectable setup (IV set, sound, occlusion)
- Wide range of delivery rate settings
- Includes drop sensor and tubing clamp
- Displays flow rate, infused volume, and total volume

**DRE Serene** »

**Portable pump compatible with most syringes**

- Built-in battery for mobile use
- LED infusion rate display
- Audible alarms



Web ID: 2105

**Also available: Refurbished Baxter syringe pumps**

**Certified Refurbished**

**Baxter Flo-Gard 6201** »

Web ID: 68

**Infusion pumps available with one or two channels**

« **Baxter Flo-Gard 6301**



Web ID: 67

**Also available: Baxter AS40, AS50, InfusO.R., and more**

**Also available: Suction Pumps**



Web ID: 1437

« **Gomco OptiVac G180**

**Portable suction unit**

- Compact and easy to use
- Power suction pump with battery backup
- Includes canister, tubing, and filter
- Additional models/sizes available

**DRE Verdana** »

**Pediatric through adult ventilation**

- High resolution display
- Integrated built-in electronic air/oxygen mixer
- Additional modes include CPAP, APRV and backup apnea
- Includes humidifier, hoses, and patient circuit
- Three levels of visual and audible alarms
- Available with optional compressor



VCV: Delivers mandatory breaths by setting tidal volume

PSV: Inspiration is usually initiated by patient effort

APRV: A time-cycled pressure mode

SIMV: Allows mandatory and spontaneous breaths

**Multiple  
Ventilation  
Modes**

« **Pulmonetic LTV-1000**

**Sophisticated,  
portable ventilator for  
use in all settings**

- Available new or refurbished
- About the size of a laptop computer
- Includes internal oxygen blender
- Integrated piston eliminates need for compressor
- Options include stand, battery pack, and display



**DRE/ORSupply.com Return Policy**  
A Customer Service Representative must authorize all returns and issue an RMA (Return Merchandise Authorization) number.

Please be prepared to provide the following information:

1. Reason for return
2. Product model number and quantity
3. Customer name and copy of sales order or invoice

- Eligibility**
- Requests for returns due to customer error will only be accepted 30 days after original shipment date.
  - A restocking fee of 40% or more will be charged. Merchandise must be properly packaged in the original shipping carton, if possible, to avoid freight damage.
  - Credit will not be given for items damaged due to improper packaging or freight damage.
  - Customers must obtain an RMA number from a DRE representative and display this number on the outside of the package. DRE reserves the right to determine in its sole discretion whether an item will be accepted as a warranty return, for quality issues, or as damaged goods.

- Charges and Handling of Credits**
- All returns are subject to inspection to determine condition of merchandise prior to issuance of credit.
  - All returns, regardless of whether the returns are determined to be legitimate, will be subject to a minimum handling/restocking fee of 40% or more of the sale price of the item. Special order items, such as upholstered chairs, are subject to higher restocking fees at the discretion of the manufacturer.
  - An additional fee of 20% of the original sale price will be charged for returned items that contain customer markings, stickers, or require repackaging to be returned to saleable condition.
  - Deposits are non-refundable.
  - Freight/shipping charges associated with the returned merchandise from the original order will not be credited.
  - Freight for returned merchandise must be shipped pre-paid freight (collect shipments will be refused).

**General Disclaimer:**  
General disclaimer in regard to battery life subject to storage conditions and temperature. Measurements and specifications provided by each manufacturer may change upon their discretion without notice. As with any piece of medical equipment, please refer to the user manual, or other manufacturer instructions, for maintenance, service requirements and details of usage. All items are subject to availability and prior sale.

**Refurbished Items:**  
Products sold used, as-is, or refurbished have no warranty unless specified, in writing, by an authorized representative of the company. Products covered by a warranty will be repaired or replaced at the discretion of the company. Purchasers should follow the Return Policy protocol above. DRE can correct any warranted defects, and may pay the cost of shipping the repaired unit back to you. Equipment that has been damaged is excluded from this warranty and will be repaired for a fee. Please contact your freight company to submit any claims for shipping damage. Please do not have your equipment repaired by a third party without prior written authorization. Doing so may void your original warranty.

All copyrights, trademarks, and brand names are the property of the company or entity to which they are listed.

Specifications subject to change without notice.

**Multiple  
Ventilation  
Modes**

Web ID: 1236

« **Puritan-Bennett 840**

**Enhanced noninvasive ventilation for neonatal through adult patients**

- Variety of software options (BiLevel, NeoMode, and more)
- Dual microprocessor and DualView touch-screen
- Supports communication with patient monitors
- Offers electronic data transfer
- Includes circuit, hoses, and mobile stand



Web ID:  
1879



**Viasys AVEA**

Non-invasive ventilation with Heliox for pediatrics and adults

Web ID:  
1880



**Viasys VELA**

Durable ventilator with pediatric through adult ventilation

Web ID:  
364

« **Siemens Servo 300A**

**Affordable ventilation for all patients**

- Auto-mode function for greater treatment possibilities
- Ventilates neonatal through adult patients
- Includes circuit, hoses, and mobile stand
- Optional graphics display available



**Dräger Babylog 8000**

**Reliable, advanced  
ventilation for neonatal  
and pediatric patients**

- Simple design and operation
- Sensitive synchronization
- Gentle but precise breathing support



## Equipment Purchasing Checklist

Following these steps can help plan your medical equipment purchase

- Before requesting a quote, make a LIST of the equipment that you need.
- If possible, include the QUANTITY of items you need.
- Indicate whether you prefer NEW or REFURBISHED equipment.
- Make a note of WHEN you will need to receive your equipment.
- Identify the equipment that you CURRENTLY USE.
- List any equipment MODELS that you prefer.
- Indicate whether you have a particular BUDGET set for your project.
- Let us know whether you would like a SHIPPING ESTIMATE.
- Please specify your facility's VOLTAGE requirements (110 or 220 volt).
- For accessories or small items, please visit [www.orsupply.com](http://www.orsupply.com)
- Make a note of your DRE salesperson's name and e-mail address.

Se habla  
Español

Find the equipment you need at

**[dremed.com/international](http://dremed.com/international)**

- Browse hundreds of new and professionally refurbished medical products
- Search by equipment type, keyword, manufacturer and product ID number
- Also shop for supplies and accessories



Check out staff and customer photos on our Facebook page!

For pricing and product information, give us a call at **502-244-6333**