

# DRE Versailles Series Accessories



# Operation Table Accessories

## Categories

Description	Page
<b>K8 Operating Table Pressure Care Mattresses</b>	<b>4</b>
Operating Table Mattresses	
<b>Head Supports</b>	<b>5</b>
Doughnut Headrings	
Ruebans Type Pillow	
Whitlock Head Support	
<b>Arm Supports</b>	<b>6</b>
Ankle/Elbow Protector	
Arm Retainers	
Carter Braine Arm Support	
Easy Arm Board	
Extended Arm Table	
Lateral Prone Arm Support	
Narrow Arm Table	
Square Arm Table	
<b>Shoulder, Spine and Torso Supports</b>	<b>8</b>
Lateral Support	
Montreal Laminectomy Support	
Pelvic Support	
Shoulder Supports	
<b>Gynecology Accessories and Supports</b>	<b>9</b>
Fluid Trap	
Knee Crutch Stirrups	
Lithotomy Pole	
Lithotomy Pole Slings	
Lloyd Davies Stirrups	
PAL Power Assisted Lithotomy Stirrups	
Perineal Instrument Tray	
<b>Leg and Foot Positioners</b>	<b>11</b>
Abduction Support	
Arthroscopic Leg Holder System	
Arthroscopic Stress Post	
D-Shaped Pillow	
The Prepper	
Total Knee Stabiliser	
Turret Heel Support	
<b>General Operating Table Accessories</b>	<b>12</b>
Abduction Wedge	
Anesthetic Screen	
Anesthetic Tube Holder	
Beanbags	
Rotary Clamps	
Rubber Ferrules	
Omni Directional Tube Holder	
Sandbags	
Stool Cushion	
Stretcher Sheets	
Handwheel and Socket Clamps	
Limb Hoist	
Foam Leg Trough	

Description	Page
Lead Hand	
Comfort Mat	
Suture Snooper	
Accessory Trolley	
<b>Pressure Relieving Gel Supports</b>	<b>17</b>
Oasis Gel Supports	
Oasis Elite Gel Supports	
Azure Gel Supports	
<b>Obstetric Pads and Mattresses</b>	<b>23</b>
Crawford Wedges	
Obstetric Mattresses	
Lithotomy Poles	
Knee Crutch Stirrups	
Stool Cushion	
<b>DRE Versailles Mobile Surgery Table Accessories</b>	<b>27</b>
Accessories - Compatible with DRE Versailles	



## Introduction

Theatre equipment specialist DRE, Inc. offers a high quality range of Operating Table Accessories. All of the items are designed and produced using the latest materials and manufacturing techniques. Our mission statement is:

‘To manufacture and market operating theatre equipment to the highest standards of materials, workmanship, innovation and after sales servicing and support; to understand our customers, their needs and moreover the solution’.

Following this policy we have developed product ranges and aspects of these products that are very unique. The products compiled in this brochure reflect our ongoing commitment to these principles.

## General Information

Pressure care: helping protect against infection  
In theatre there is a strong emphasis on pressure care – minimising a patient’s risk of developing pressure sores caused by immobility during an operation and the use of anesthesia. Pressure sores are of particular concern in tissue viability as broken skin inevitably makes the patient more vulnerable to infection.



We have products specifically designed for pressure care:

K8 mattresses – available for most makes of operating table and fitted as standard on DRE Versailles mobile surgery tables:

- made of a gel-feel foam which moulds to the patient’s body (in the science of Kinetics, there are eight points of contact the human body makes when lying on a hard surface – hence K8), thus relieving the pressure by distributing the patient’s weight as evenly as possible
- antibacterial and microbial fungistatic properties - the outer cover of each mattress is made of an anti-microbial, water resistant, vapour permeable material which allows the skin to breathe. The unique fabric sealing system effectively welds the outer material together to give an almost unbreakable seal, making the K8 a class leader in infection control.

Azure and Oasis gel pad and positioning products:

- Azure uses two unique silicone materials which replicate skin and body tissue – it conforms to individual’s body contours to provide maximum pressure relief, shear reduction and protection against nerve damage
- Oasis and Oasis Elite ranges consist of both silicone gel pads and positioning devices, which evenly distribute body weight and don’t ‘bottom out’ – qualities essential for the relief of pressure during surgery

## Guaranty

All of our products are manufactured to documented in-house quality control procedures, individually inspected and guaranteed against any defects in material or workmanship.

## After Sales Service


Our qualified sales team is available for consultation and can often provide samples on site for appraisal.



# K8 Operating Table Pressure Care Mattresses

## K8 Pressure Care Mattresses


Recent evidence suggests that any patient undergoing surgery for more than three hours faces a risk of developing pressure sores, caused by the immobility during the operation and the use of anesthesia. These products are designed to minimize these effects: the K8 Pressure Care Mattress is available for all makes of operating table and can be custom made to specific requirements.



The mattresses are made of a gel-feel foam which moulds to the patient's body, in the science of Kinetics, there are eight points of contact the human body makes when lying on a hard surface – hence K8, thus relieving the pressure by distributing the patient's weight as evenly as possible and minimising the risk of decubitus ulcers.



- High Risk Pressure Care Mattress
- Independently Tested
- Latex Free
- Antibacterial/ Vapour Permeable
- Welded Seams
- X-Ray Translucent
- Specialist Fitting Service Available




Ref No.	Description
10005	M Type Table Mattress K8 Pressure Care
10010	MR Type Table Mattress K8 Pressure Care
10015	J3 Type Table Mattress K8 Pressure Care
10020	RX500 Type Table Mattress K8 Pressure Care
10025	RX600 Type Table Mattress K8 Pressure Care
10030	MXP Type Table Mattress K8 Pressure Care
10050	Belfast Type Table Mattress K8 Pressure Care



# K8 Operating Table Pressure Care Mattresses




Ref No.	Description
10070	ALM Type Table Mattress K8 Pressure Care
10075	Maquet Alphastar Type Table Mattress K8 Pressure Care
10077	Maquet Betastar Type Table Mattress K8 Pressure Care
10078	Maquet Orthostar Type Table Mattress K8 Pressure Care
10080	T20 Type Table Mattress K8 Pressure Care
10085	Steris Type Table Mattress K8 Pressure Care



All mattresses are supplied with polyurethane Velcro fastenings as standard. Specialist AI-Grips are available on request.

Mattresses can be supplied to your own dimensions. Our qualified sales team are available for on-site consultations for specific requirements.

## Head Supports



### Head Supports

All specialisms including ENT, dental, ophthalmics, facial maxillary and neurosurgery are covered by a range which includes head supports for patients in prone and supine positions.



Ref No.	Description
10170	Doughnut Heading Adult



Ref No.	Description
10175	Doughnut Heading Child



Ref No.	Description
10180	Ruebans Pillow Adult



Ref No.	Description
10185	Whitlock Horseshoe Headrest






## Arm Supports

### Arm Supports

We offers a wide selection of arm supports including hand tables, for positioning and protecting patients in prone, supine or lateral positions.



Ref No.	Description
10105	Arm Retainer L Shaped
10110	Arm Retainer J Shaped
10112	Arm Retainer Extra Wide

Ref No.	Description
10115	Arm Retainer Pad Suitable for 10105 & 10110
10117	Arm Retainer Pad Extra Wide

Ref No.	Description
10152	Elbow Board
10153	Elbow Board Pad



Ref No.	Description
10195	Carter Braine Arm Support Complete Requires clamp 10250 (x1)
10196	Carter Braine Arm Support Pad Only
10198	Carter Braine Arm Support 100mm Pan Complete Requires clamp 10250 (x1)
10199	Carter Braine Arm Support Pad Only 100mm Pan



# Arm Supports



- | Ref No. | Description  |
|---------|--|
| 10220   | Narrow Arm Table Complete 590mm x 210mm<br>Requires Clamp 10250 (x1) |
| 10221   | Narrow Arm Table Pad Only 590mm x 210mm                              |



- | Ref No. | Description  |
|---------|--|
| 10230   | Extended Arm Table Complete (X-Ray Translucent)<br>Requires Clamp 10250 (x2) |
| 10231   | Extended Arm Table Pad Only  |




- | Ref No. | Description                                 |
|---------|---|
| 10774   | Easy Arm Board Lift-to-Position (Euro Rail) |
| 10775   | Easy Arm Board Lift-to-Position (UK Rail)   |



- | Ref No. | Description                             |
|---------|---|
| 10225   | Square Arm Table Complete 450mm x 415mm |
| 10226   | Square Arm Table Pad Only 450mm x 415mm |



- | Ref No. | Description               |
|---------|---------------------------|
| 10773   | Lateral/Prone Arm Support |



# Shoulder, Spine and Torso Supports

## Spine and Torso Supports

An extensive range of shoulder, spine and torso supports with a combined emphasis on patient safety and surgical access. It has recently been complemented by the shoulder arthroscopy attachment for the DRE Versailles Mobile Surgery Table.



Ref No. Description  
10160 Montreal Laminectomy Support



Ref No. Description  
10787 Shoulder Supports  
Secures patient in trendelenberg



10789 Weightless Shoulder Suspension System  
10791 Disposable Arm Traps (Case of 6) Sterile



Ref No. Description  
10200 Lateral Support Complete  
Requires Clamp 10250 (x1)  
10201 Lateral Support Pad  
10205 Pelvic Support Complete  
Requires Clamp 10250 (x1)  
10206 Pelvic Support Pad Only




Ref No. Description  
21340 Shoulder Arthroscopy Support  
For use with Versailles Trolley System



Ref No. Description  
10770 Lateral Braces, 3 piece  
Individual pieces available on request

# Gynecology Accessories and Supports



## Gynecology Accessories and Supports

From the powered Pal Pro™ stirrups to conventional lithotomy poles, DRE offers every accessory for this surgical speciality.



Ref No.	Description
10210	Lithotomy Poles Complete Requires Clamps 10245 or 10246 (x2)
10211	Lithotomy Poles Sleeves Only
10212	Lithotomy Poles Padded Slings Only
10213	Lithotomy Poles Canvas Slings Only



Ref No.	Description
10235	Perineal Instrument Tray
10236	Multi-pin for Maquet
10237	Multi-pin for Eschmann
10238	Multi-pin for Schearer
10239	Multi-pin for Steris




Ref No.	Description
10215	Lloyd Davies Stirrups Complete Requires Clamps 10245 or 10246 (x2)
10216	Lloyd Davies Stirrups Pads Only



Ref No.	Description
10765	Knee Crutch Stirrups
10767	Fluid Trap - Euro Catcher System
10768	Drape, Disposable Sterile (pack of 10)
10769	Drape, Disposable Non Sterile (pack of 10)





## Gynecology Accessories and Supports



### Pal Pro™ Stirrups

The Pal Pro™ stirrup is our latest in power-assisted legholders. This stirrup allows control of the level of lithotomy and abduction from the release handle. Because of the unique piston design, the stirrup does most of the work lifting the patient's legs. The floating boot is self adjustable, to minimise pressure on the calf. Positioning the patient has never been easier.

- Unhindered surgical access throughout full lithotomy range
- Abduction adjustments made quickly from end of boot
- Floating boot design protects superficial nerves
- FootHugger™ boot pads provide comfortable and secure leg positioning
- Easy-to-squeeze handle
- Intra-operative repositioning by one person from end of boot
- Simplified setup reduces staffing requirements and time
- Feather Lift power assist minimises staff lifting
- Easy-to-clean boot and protective piston cover
- Accommodates patient weight up to 159 kg

### Heavy Duty Stirrup


A new addition to our range is the Heavy Duty Stirrup which easily accommodates patients with larger calves.

The reinforced rod and joint assembly and the flexible boot make the system capable of positioning patients weighing up to 363 kg

### Mobile Storage Cart

The perfect way to store and protect your investment, this strong and lightweight cart holds and transports the stirrups or legholders.

Compact and state-of-the-art, it can hold any stirrup from any manufacturer. The sturdy moulded four-castor polymer base is dent-resistant and the deep storage bin is a great place to store clamps and pads securely.



Also Available:  
Heavy Duty Stirrup capable of positioning patients weighing up to 363 kg



Ref No.	Description
10750	Pal Pro Power Assisted Stirrups
10758	Power Assisted Heavy Duty Stirrups
10760	Handy Kart

#### Accessories:

Ref No.	Description
10710	FootHugger Boot Pads for Pal Pro Supports
10740	Boot Pads Gel for Pal Pro Supports
10761	Handwheel Clamp, Universal - Single Unit




# Leg and Foot Positioners



## Leg and Foot Positioners

We can offer a wide range of leg and foot positioners from arthroscopic leg holders and total knee stabilisers to abduction board and heel supports.



- | Ref No. | Description                          |
|---------|--------------------------------------|
| 10145   | D Shaped Pillow 480mm x 125mm x 65mm |
| 10150   | Turret Heel Support Single           |



- | Ref No. | Description                    |
|---------|--------------------------------|
| 10780   | Arthroscopic Leg Holder System |



- | Ref No. | Description                              |
|---------|--|
| 10155   | Abduction Support for 75mm Deep mattress |



- | Ref No. | Description                    |
|---------|--------------------------------|
| 10785   | Total Knee Stabiliser Complete |



- | Ref No. | Description             |
|---------|-------------------------|
| 10776   | Athroscopic Stress Post |



- | Ref No. | Description |
|---------|-------------|
| 10798   | The Prepper |

# General Operating Table Accessories

## General Accessories

Specialisms aside, there are huge number of general operating table accessories to be found in DRE's portfolio – from tube holders and screens to clamps, beanbags and stool cushions.



Ref No.	Description
10120	Anesthetic Tube Holder Flat
10122	Anesthetic Tube Holder Ang.



Ref No.	Description
10190	Anesthetic Screen Complete Requires Clamp 10250 (x1)
10191	Anesthetic Screen Only Requires Clamp 10250 (x1)
10192	Anesthetic Screen Sleeve Only



Ref No.	Description
10125	Omnidirectional Tube Holder Functions as flexible tube holder & Anesthetic Screen



Ref No.	Description
10130	Stool Cushion 330mm dia. Stool Cushions can be supplied to your own dimensions.

# General Operating Table Accessories




Ref No. Description  
 10240 Limb Hoist  
 Assists user during limb preparation




Ref No. Description  
 10248 Rotary Clamp,  
 Universal  
 Suitable for BS or  
 Euro Rails



Ref No. Description  
 10250 Handwheel Socket  
 Clamp, Universal  
 Suitable for BS or Euro Rails



Ref No. Description  
 10245 Easy Lock Rotary  
 Clamps (UK Rail)  
 10246 Easy Lock Rotary  
 Clamps (Euro Rail)



Ref No. Description  
 10761 Handwheel Clamp,  
 Universal  
 Suitable for BS or  
 Euro Rails



Simply snap on the clamp to any rail location. Accepts any round post up to 18mm. Quick and easy pull-to-release lever. Lithotomy wheel for accessory angle adjustment.



Ref No. Description  
 10763 Easy Lock Blade Clamp  
 (UK Rail)  
 10764 Easy Lock Blade Clamp  
 (Euro Rail)



# General Operating Table Accessories




- | Ref No. | Description                  |
|---------|------------------------------|
| 10255   | Sandbag Small 200mm x 150mm  |
| 10256   | Sandbag Medium 300mm x 200mm |
| 10257   | Sandbag Large 380mm x 250mm  |
- Sandbags can be supplied to your own dimensions.



- | Ref No. | Description  |
|---------|--------------|
| 10796   | Suture Snoop |



- | Ref No. | Description           |
|---------|-----------------------|
| 10265   | Ankle/Elbow Protector |



- | Ref No. | Description |
|---------|-------------|
| 10797   | Comfort Mat |



- | Ref No. | Description                       |
|---------|-----------------------------------|
| 10270   | Rubber Ferrules<br>(pk of 5) 19mm |
| 10271   | Rubber Ferrules<br>(pk of 5) 25mm |





- | Ref No.  | Description              |
|----------|--------------------------|
| L6410-37 | Lead Hand Splint - Adult |

# General Operating Table Accessories



Ref No.	Description
10275	Beanbag Vacuum Full Body 1000mm x 730mm
10276	Beanbag Vacuum Square 500mm x 400mm
10277	Beanbag Vacuum Rectangular 460mm x 230mm
10278	Beanbag Vacuum Head Stabiliser 360mm x 230mm
10279	Beanbag Vacuum Full Length 1650mm x 720mm



Ref No.	Description
10285	Foam Leg Trough (covered)
10290	Abduction Wedge



10550	Crawford Wedge Small 550mm x 300mm x 100mm
10555	Crawford Wedge Large 550mm x 510mm x 200mm

# General Operating Table Accessories

## Accessory Trolley



This ergonomically designed stand offers a complete storage solution for general operating theater and Versailles accessories.

Adaptable, robust and multi-functional, it is ideal for modern theatre environments and clinical areas.

Specific storage for Versailles accessories with locations for:

- 4 cotsides
- 2 leg sections
- power cable
- 2 transfusion poles

The trolley still provides positions for a large number of general accessories.




Ref No.	Description
21370	Accessory Trolley British Standard Rail
21370E	Accessory Trolley European Standard Rail
21371	Accessory Trolley Spare Trough



### General Specifications: Versailles Accessory Stand

Frame	Teflon <sup>®</sup> Coated		
Base	Acrylic Capped, Vacuum Formed Base Antibacterial		
Castors	100mm Twin Wheel Polyurethane Tyred Non-Corrosive Castors		
Dimensions	Height	Width	Depth
Stand	1370mm	720mm	720mm
Storage Box	280mm	680mm	230mm
Max. Load	100kg		

# Pressure Relieving Gel Supports



## Pressure Relieving Gel Supports



# Pressure Relieving Gel Supports

## Pressure Relieving Gel Supports

Oasis offers a complete range of lightweight, reusable and superior pressure relieving pads and positioners for use in the operating theatre. The range has been designed to re-distribute weight and provide excellent pressure relief at all risk levels.

It is a fact that up to 25% of bed sores begin in surgery and so it is important to ensure that every operating theatre is equipped with the most effective pressure relieving products.



### OASIS:

- Silicone superior pressure relief
- Conforms to body shape
- Latex free and skin compatible
- Comprehensive range
- Can be heated or cooled



### OASIS ELITE:

- Foam and silicone gel combination
- Enhanced protection
- Lightweight
- Highly durable
- Easy to clean and maintain



### AZURE:


- Revolutionary new gel combination
- Contoured protection
- Redefining pressure care
- Advanced surgical positioning and protection
- 12 month guarantee

## Oasis<sup>®</sup> Gel Supports

The Oasis range of reusable silicone gel pads and positioners provide superior pressure relief to the patient during various surgical procedures.

















The gel evenly distributes body weight, essential for the relief of pressure during surgery.

All Oasis products can be heated or cooled, have an X-ray attenuation of 89.7%, are latex free and non-allergenic. They conform to the patient's body shape and have a 12 month guarantee against manufacturing defects.



Ref No.	Description	
10801	Closed Head Ring, Pediatric 90mm x 40mm x 20mm	
10802	Closed Head Ring, Adolescent, 140mm x 60mm x 35mm	
10803	Closed Head Ring, Adult, 200mm x 75mm x 45mm	
10804	Open Head Ring, Pediatric 90mm x 40mm x 20mm	
10805	Open Head Ring, Adolescent 140mm x 60mm x 35mm	
10806	Open Head Ring, Adult 200mm x 75mm x 45mm	
10807	Operating Table Pad 520mm x 520mm x 10mm	
10808	Operating Table Pad 1150mm x 520mm x 10mm	
10809	Operating Table Pad 1800mm x 520mm x 10mm	
10811	Hip Pad 520mm x 520mm x 10mm	

Ref No.	Description	
10812	Table Pad, Head Section 500mm x 225mm x 10mm	
10813	Table Pad, Torso Section w/ Per. cut-out 500mm x 990mm x 10mm	
10814	Table Pad, Foot Section 500mm x 500mm x 10mm	
10815	Table Pad Set, 3 Piece Consists of 1 each - OA036, OA037, OA038	
10828	Ulnar/Brachial Nerve Protector 380 x 115 x 8mm	
10829	Ulnar/Brachial Nerve Protector 460 x 210 x 8mm	
10830	Armboard Pad 390mm x 115mm x 10mm	
10831	Large Armboard Pad 600mm x 200mm x 10mm	
10832	Armboard Pad, Long 600mm x 125mm x 10mm	
10833	Contoured Armboard Pad, 510mm x 150mm x 35mm	

Ref No.	Description		Ref No.	Description	
10816	Perineal Table Pad 520mm x 520mm x 10mm		10834	Arm Retainer Cover 170mm x 110mm x 10mm	
10817	Profile Pad, 520mm x 520mm x 12mm		10835	Arm Retainer Cover, Wide 170mm x 132mm x 10mm	
10818	3/4 Length Profile Pad 1145 x 520 x 12mm		10836	Ankle Strap Stirrups x 2 310mm x 45mm x 10mm	
10819	3/4 Length Profile Pad w/ cutout, 1145mm x 520mm x 12mm		10837	Large Ankle Strap Stirrups 440mm x 45mm x 10mm	
10823	Patient Positioner, Ex-Small, 160mm x 60mm x 45mm		10838	Crutch Stirrup Pads x 2 315mm x 235mm x 10mm	
10824	Patient Positioner, Small, 330mm x 70mm x 50mm		10839	Crutch Stirrup Pads Large x 2 395mm x 245mm x 10mm	
10825	Patient Positioner, Medium, 300mm x 120mm x 60mm		10840	Crutch Stirrup Pads Large x 2 520mm x 300mm x 10mm	
10826	Patient Positioner, Large, 360mm x 125mm x 70mm		10851	Perennial Post Pad 300mm x 260mm x 10mm	
10841	Lithotomy Pole Sleeves x 2		10852	Traction Boot Pad 375mm x 245mm x 5mm	
10842	Patient Positioning Belt 505mm x 85mm x 8mm		10853	Heel and Ankle Protector 230mm x 200mm x 20mm	
10845	Heel Pads x 2 200mm x 100mm x 65mm		10854	Chest Roll w/ Flat Bottom 305mm x 100mm x 75mm	
10846	Heel Cups x 2		10855	Chest Roll w/ Flat Bottom 405mm x 100mm x 75mm	
10848	Heel Cups, Case of 32 Pairs		10856	Chest Roll w/ Flat Bottom 510mm x 100mm x 75mm	
10849	Face Pad 165mm x 80mm x 8mm		10858	Oasis First Aid Kit	

# Pressure Relieving Gel Supports




## Oasis Elite

Oasis Elite products combine a unique lightweight foam with silicone Trugel™ to provide superior pressure relief and enhanced patient protection in theatre.

The range conforms to the patient's body shape, distributing weight over a greater surface area to reduce potential nerve damage and has undergone rigorous durability tests.

Oasis Elite products are more than 50% lighter than comparable gel-only devices. The products are easy to clean and maintain.



Ref No.	Description	Image	Ref No.	Description	Image
10861	Closed Head Ring, Adolescent, 140mm x 34mm		10873	Supine Head Rest, Adult 272mm x 237mm x 90mm	
10862	Closed Head Ring, Adult, 200mm x 50mm		10874	Supine Head Rest, Adolescent 270mm x 230mm x 60mm	
10863	Open Head Ring, Adolescent 140mm x 134mm		10875	Heel Pads x 2 180mm x 100mm x 70mm	
10864	Open Head Ring, Adult, 200mm x 54mm		10876	Dome Positioner 350mm x 150mm x 75mm	
10865	Table Pad, 520 x 520 x 20mm		10877	Chest Roll, 305mm x 150mm x 100mm	
10866	3/4 Table Pad 1150mm x 520 x 20mm		10878	Chest Roll, 405mm x 150mm x 100mm	
10867	Armboard Pad, 400mm x 125mm x 20mm		10879	Chest Roll, 510mm x 150mm x 100mm	
10868	Armboard Pad 600mm x 125mm x 20mm		10880	Chest Roll, 405mm x 150mm x 150mm	
10869	Contoured Armboard Pad, 510mm x 150mm x 35mm		10881	Chest Roll, 510mm x 150mm x 150mm	
10870	Contoured Armboard Pad 400mm x 150mm x 35mm		10882	Oasis Elite Lateral Positioner 720mm x 500mm x 160mm	
10871	Prone Head Rest, Adult 280mm x 240mm x 140mm		10883	Elite Prone Chest Positioner 600mm x 120mm x 520mm	
10872	Prone Head Rest, Adolescent 280mm x 240mm x 130mm				




# Pressure Relieving Gel Supports

## Azure® Gel Supports

Azure has come about as a result of Trulife's commitment to continually provide improved levels of pressure care products to both patients and healthcare systems.

The revolutionary new gel combination works together to replicate and act as an additional layer of protective tissue and skin, conforming to the individual user's shape in a more contoured and less structured manner. This results in effective protection against nerve damage as well as pressure sore prevention.

Azure products are re-defining pressure care and are set to change the standards in theatre.



Ref No. Description  
10885 Head Protector  
260mm x 250mm x 50mm



10886 Profile Pad,  
520mm x 520mm  
x 10mm



10887 Profile Pad,  
1150mm x 520mm x 10mm



10888 Sacral Protector  
400mm x 380mm x 25mm

10889 Arm Protector, Short  
330mm x 140mm x 30mm




10890 Arm Protector,  
Long  
460mm x 180mm x 35mm




10891 Heel Protectors x 2  
190mm x 175mm x 30mm



10892 Leg Protector  
515mm x 180mm x 40mm



10893 Patient Positioner  
410mm x 145mm x 60mm



## Obstetric Pads and Mattresses



## Obstetric Pads and Mattresses

# Obstetric Pads and Mattresses

## Obstetric Pads and Mattresses

The K8 mattresses manufactured from gel feel foam actually moulds to the patient's body, thus relieving pressure by distributing the patient's weight as evenly as possible. Antibacterial, water resistant materials, which are vapour permeable and latex free, allow the skin to breathe. The unique fabric sealing system effectively welds the outer material together to give an almost unbreakable seal making the K8 a class leader in infection control.




Ref No.	Description
10500	Birthday Bed Type Complete Delivery Bed Mattress
10515	Birtright Type Complete Delivery Bed Mattress
10520	Lic Optima Type Complete Delivery Bed Mattress
10525	Hoskins Steel Type Complete Delivery Bed Mattress
10530	EVAB Type Delivery Bed Mattress

Mattresses can be supplied to your own dimensions. Our qualified Sales Team are available for on-site consultation for any specific requirements.



Ref No.	Description
10210	Lithotomy Poles Complete Requires Clamp 10245 or 10246 (x2)
10211	Lithotomy Poles Sleeves Only
10212	Lithotomy Poles Padded Slings Only
10213	Lithotomy Poles Canvas Slings Only
10765	Knee Crutch Stirrups



Ref No.	Description
10130	Stool Cushion 330mm dia.
10550	Crawford Wedge Small 550mm x 300mm x 100mm
10555	Crawford Wedge Large 550mm x 510mm x 200mm

# DRE Versailles Series Mobile Surgery Table Accessories



## Accessories for DRE Versailles Series Mobile Surgery Tables

# DRE Versailles Mobile Surgery Table Accessories

## DRE Versailles Mobile Surgery Table

With both powered and manual versions available, this combination of patient transport trolley and operating table is a complete multi-purpose mobile operating system. The DRE Versailles Mobile Surgery Table is a familiar sight in operating theatres, day treatment centers, and increasingly in private clinics and GP surgeries.



Ref No.	Description
21300	DRE Versailles P100 — Powered Functions
21310	DRE Versailles M75 — Manual Functions

### Powered Functions:

- Electrically powered functions
- Superior Surgical Access
- Traversing Patient Platform
- 100% 'C' Arm Access
- Lateral Tilt
- 250kg Weight Limit

### Manual Functions:

- Mechanically Powered Functions
- Increased Surgical Access
- Supremely Lightweight and Maneuverable
- 75% 'C' Arm Access
- Low Cost Solution
- 160kg Weight Limit



Ref No.	Description
10190	Anesthetic Screen Complete
10191	Anesthetic Screen Only Requires Clamp 10250 (x1)



Ref No.	Description
10195	Carter Braine Arm Support Pad Complete Requires Clamp 10250 (x1)



Ref No.	Description
10200	Lateral Support Complete Requires Clamp 10250 (x1)



Ref No.	Description
10205	Pelvic Support Complete Requires Clamp 10250 (x1)



Ref No.	Description
10210	Lithotomy Poles Complete Requires Clamp 10245 or 10245 (x2)

# DRE Versailles Series Mobile Surgery Table Accessories



Ref No. 10220  
 Description  
 Narrow Arm  
 Table Complete  
 590mm x 210mm  
 Requires Clamp  
 10250 (x1)



Ref No. 10750  
 Description  
 Pal Pro Power Assisted Stirrups



Ref No. 10225  
 Description  
 Square Arm  
 Table Complete  
 450mm x 415mm  
 Requires Clamp  
 10250 (x2)



Ref No. 21280  
 Description  
 Emergency  
 Venturi Suction  
 System



Ref No. 10230  
 Description  
 Extended Arm  
 Table Complete  
 (X-Ray  
 Translucent)  
 Requires Clamp  
 10250 (x2)



Ref No. 21320  
 Description  
 Ophthalmic  
 Head Support  
 c/w Sculptured Pad



Ref No. 10245  
 Description  
 Easy Lock  
 Rotary Clamps  
 (B/S Rail)



Ref No. 21322  
 Description  
 Full Width  
 Articulating  
 Head Rest

Ref No. 10246  
 Description  
 Easy Lock  
 Rotary Clamps  
 (Euro Rail)



Ref No. 10250  
 Description  
 Handwheel  
 Socket Clamp,  
 Universal



Ref No. 21325  
 Description  
 Full Width  
 Non Articulating  
 Head Rest

# DRE Versailles Mobile Surgery Table Accessories



Ref No. Description  
21340 Shoulder Arthroscopy Support



Ref No. Description  
21358 Perineal Instrument Tray

Ref No. Description  
21359 Abduction Support For Varicose Vein Surgery






Ref No. Description  
21365 Oxygen Delivery System  
Requires Clamp 10250 (x1)



Ref No. Description  
21350 Padded Cotside Covers, Grey  
21351 Padded Cotside Covers, Child Print

Ref No. Description  
21360 K8 Pressure Care Mattress



Ref No. Description  
21355 Articulating Removeable Leg Section  
21357 Removeable Foot Section c/w Full Length Side Rail



Ref No. Description  
21390 Push/Pull Bar  
21395 Foot End Extension For Surgery System



# DRE Versailles Mobile Surgery Table Accessories

This ergonomically designed stand offers a complete storage solution for general operating theatre and Versailles accessories.

Adaptable, robust and multi-functional, it is ideal for modern theatre environments and clinical areas.



Ref No.	Description
21370	Accessory Trolley British Standard Rail
21370E	Accessory Trolley European Standard Rail
21371	Accessory Trolley Spare Trough

## Delivering excellence in customer service

We are committed to providing the highest possible standards of after care service and product reliability. The product range includes:

- Versailles Surgery systems
- Quality Theatre Furniture
- Tourniquet equipment
- Surgical Instruments
- Electro Surgical Accessories



# DRE Versailles Series Accessories

