

# All-in-one ENT System

with Futura Evo Power Procedure Chair

The Euroclinic Otocompact Luxe is an all-in-one ENT procedure system. Its modular construction and mobility replace built-in cabinetry and inefficiently designed work spaces. The Luxe combines all of the commonly used elements of an outpatient setting with cutting-edge technology for superior performance.

The Luxe's Italian design provides an excellent experience for both the patient and the doctor. With a top of the line ENT chair and precision-made mobile cabinet housing constructed from the highest quality materials, the Luxe is ergonomic, attractive, and provides a sophisticated look that fits well in large or small spaces.

The Luxe integrates wireless controls, instrument table, easily-maneuverable microscope, LED light source, and endoscopy video system into one easy-to-use unit. The Luxe is the ideal option for new clinics or doctors replacing several items. It is unique to the market, high performing, and very cost effective. It is a truly complete office solution that is engineered to significantly increase productivity.



**Zeiss Opmi Pico  
Microscope**



**Futura Evo  
Procedure Chair**



**LED Light Source**



# Otocompact Professional Luxe

## ENT Workstation

Equipment for the way *you* operate

### Key Components

#### Integrated Italian Leather Seating

- » Futura Evo ENT chair with four powerful motors, three electronic memory settings, wireless foot control, and non-deformable foam padding
- » SPD/A doctor's stool built from stainless steel with height adjustment and swiveling back rest

#### Mobile European Styled Workstation

- » Instrument management system with electronic console for dual control of chair and instruments
- » Powerful suction system with dry pump
- » Insufflation system with nasal olive and autoclavable cannula
- » Auricular washing syringe with anti-sprinkle shield and microprocessor-controlled heating element
- » Larynx mirror heater programmable for 10-120 seconds
- » Air compression system for medications with 2 atomizers for liquids and powders
- » Antimicrobial filter is 99.99% effective against viruses and bacteria

#### German Engineered ENT Microscope

- » Zeiss Opmi-PICO otomicroscope with precision mount and endoport for camera connection

#### Advanced Lighting System

- » Bright LED light source (180 Lumen at 5,500° K)
- » Neutral LED light source (180 Lumen at 5,000° K) with Storz connector

#### Integrated Video Endoscopy Suite

- » Wireless Camera EVS ED410 with ergonomic design and wireless foot pedal for hands-free photos and video
- » 23" HP TouchSmart PC monitor on articulating arm with wireless keyboard and mouse
- » Video capture and data management software
- » Reader/Writer smart card kit for easy record keeping



### Futura Evo Specifications

Power Supply Tension	110V or 230V-50/60Hz
Linak Motor Low Tension	24 Vdc
Motor Lifting Push	6000N fast run
Back Rest Motor Push	6000N
Seat Motor Push	6000N
Unique Leg Rest Motor Push	6000N
Lifting Capacity	180 kg
Absorbed Power	120 VA
Weight	120 kg
Base Length	800 mm
Maximum Height of Leg Rest	1250 mm

Specifications subject to change without notice.



### Otocompact Professional Evo Specifications

Power Supply Tension	110/120V or 220/240 V
Frequency	50/60Hz
Absorbed Power (MAX configuration)	1200 VA
Weight	120 kg
Maximum Number of Instruments	4
Maximum Air Capacity	26/min
Dry Suction Pump	-86 Kpa (645 mm Hg, 50 lt/min)
Oil Suction Pump	-86 Kpa (645 mm Hg, 70 lt/min)
Power Supply Fuses	F1, F2=6, 3 AT
Fuse External Outlet	F3=4 AT